

CelGel 1kb DNA ready ladder (100-10 000bp)

CM1340-0500 100 gel lanes

Store the kit at -20°C



CELTIC MOLECULAR

Explore . Discover . Quality

What's in the box

Cat no.	CM1340-0500
Pack size	100 gel lanes
CelGel 1kb DNA ready ladder	500µl

Product description

A molecular weight marker designed to determine the weight and size of double stranded DNA during gel electrophoresis analysis. The CelGel 1kb DNA ready ladder contains 13 unique fragments ranging from 100 to 10 000bp in size, consisting of PCR products and restriction enzyme digested plasmids. The approximate mass of each DNA band (when loading 5µl per well) is supplied in below annotated gel image. This enables the estimation of the mass of sample DNA present in bands of similar intensity and size. Note that the 1kb and 3kb fragments have increased intensity to serve as easy identifiable reference points.

Shipping and storage instructions

Store the kit at -20°C for up to 24 months. This product is shipped on ice blocks and can be kept at 4°C for 12 months without affecting the product performance.

Notes

For research use only.

Technical support

For technical support please e-mail info@celticmolecular.com

Properties

Proprietary plasmids digested with appropriate restriction enzymes are combined with PCR amplicons and purified via phenol extraction. Supplied in a solution 10mM Tris-HCl (pH 8.0), 10mM EDTA and combined with loading dye.

Size range : 100 – 10 000bp
Number of bands : 13
Concentration: 50µg/500µl
Tracking dyes: Xylene Cyanol FF and Bromophenol Blue
Format: Ready to load

Reaction set-up

Load 5µl ready ladder per well (contains loading dyes Xylene Cyanol FF and Bromophenol Blue).

Tip: use with 1µl CelGel ready DNA stain and loading dye (CM1410-1000) for instant visualization of DNA bands on a UV or visible -light transilluminator.

Visual weight and size distribution

